

**CAUTION: UNITED STATES FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN TRAINED IN THE USE OF THE CERENE CRYOTHERAPY DEVICE**

## INDICATIONS FOR USE

The *Cerene Exhaust Collection Bag (Exhaust Collection Bag)* is indicated for use with the *Cerene Cryotherapy Device*. The *Exhaust Collection Bag* collects the nitrous oxide (N<sub>2</sub>O) gas exhausted from the *Cerene Cryotherapy Device* during the endometrial cryoablation procedure and facilitates disposal of the gas when the procedure is completed.

## DEVICE DESCRIPTION

The *Exhaust Collection Bag* is made of a polyurethane blend and is provided non-sterile and individually labeled in boxes of ten. The *Exhaust Collection Bag* has a port that should be connected with the *Exhaust Hose* of the *Cerene Cryotherapy Device (Cerene Device)* prior to the procedure. The *Exhaust Collection Bag* can be emptied via wall suction once the procedure is complete by using the wall suction adapter provided in the box of ten *Exhaust Collection Bags*. The *Exhaust Collection Bag* can also be emptied outdoors by opening the plug at the bottom of the bag and rolling the bag down to release the exhaust gas.

## WARNINGS

- The *Exhaust Hose* must not be occluded during treatment. *Exhaust Hose* occlusion will cause the *Cerene Device* to shut down to avoid excessive uterine pressures and may result in early termination of the procedure.
- Always connect the *Exhaust Hose* to the *Exhaust Collection Bag* prior to turning on the *Cerene Device*. Connecting the *Exhaust Hose* to the *Exhaust Collection Bag* during the procedure could result in a temporary occlusion, could cause the *Cerene Device* to shut down to avoid excessive uterine pressures, and may result in early termination of the procedure.
- Modification of the *Exhaust Collection Bag* may result in serious patient injury.

## PRECAUTIONS

- Inhalation of nitrous oxide gas (exhaust gas) can cause short-term effects such as breathing difficulty, drowsiness, headache, and asphyxia. Use of the *Exhaust Collection Bag* accessory is required.
- Chronic exposure to nitrous oxide gas (exhaust gas) can lead to adverse reproductive, neurological, and hematological effects. Use of the *Exhaust Collection Bag* accessory is required.

## PERFORMANCE CHARACTERISTICS

The *Exhaust Collection Bag* provides a means of containing and safely disposing of the nitrous oxide gas from a *Cerene* procedure in order to minimize the occupational risks associated with long-term, nitrous oxide exposure.

## SERIOUS INCIDENT REPORTING

In the European Union, any serious incident that has occurred in relation to the device should be reported to Channel Medsystems (+1-510-338-9301 or [support@cerene.com](mailto:support@cerene.com)) and the competent authority of the Member State in which the user is established. To report a product complaint or adverse event in the United States, please contact Channel Medsystems (+1-510-338-9301 or [support@cerene.com](mailto:support@cerene.com)) or FDA at 1-800-FDA-1088 or [www.fda.gov/medwatch](http://www.fda.gov/medwatch)



**Figure 1. Exhaust Collection Bag (an accessory to the Cerene Cryotherapy Device), shown here hanging from an IV pole**

## GENERAL INSTRUCTIONS

- Hang the *Exhaust Collection Bag* on an IV pole.
- Before the *Cerene Device* is turned on, connect the end of the *Exhaust Hose* to the *Exhaust Collection Bag* until you hear an audible click.

**NOTE:** Ensure that the *Exhaust Hose* is not kinked or obstructed at any time.

**NOTE:** Do not connect the *Exhaust Hose* to a scavenging system or to wall suction directly.

## DISPOSAL

- Once the *LCD Screen* indicates the *Cerene Device* is safe for disposal, disconnect the *Exhaust Collection Bag* from the *Cerene Device*, and dispose of the *Cerene Device* per local procedure.
- Empty the *Exhaust Collection Bag* and dispose per custom and practice.
  - Outdoors: Open the plug, hold the bag below waist level, and roll the bag from the top down.
  - Via wall suction: Attach the wall suction adapter found in the box of *Exhaust Collection Bags* to the mating fitting at the bottom of the bag. Attach the adapter to wall suction.

## PACKAGING AND STORAGE

The *Exhaust Collection Bag* is provided non-sterile and individually packaged, or individually packaged in boxes of ten. Do not use the *Exhaust Collection Bag* if the package appears damaged in any way. Do not reuse the *Exhaust Collection Bag*. Re-use of the *Exhaust Collection Bag* could promote cross-contamination. Store the *Exhaust Collection Bag* at -18°C – 55°C.

## SYMBOL DEFINITIONS

Symbol	Standard/Symbol Reference Number	Title of Symbol	Definition
	ISO 15223-1, Clause 5.1.6	Catalog number	Indicates the manufacturer's catalog number so that the medical device can be identified.
	ISO 15223-1, Clause 5.1.1	Manufacturer	Indicates the medical device manufacturer.
	ISO 15223-1, Clause 5.4.3	Consult Instructions for Use	Indicates the need for the user to consult the instructions for use.
	ISO 15223-1, Clause 5.4.2	Do not reuse	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	ISO 15223-1, Clause 5.3.7	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.
	ISO 15223-1, Clause 5.1.5	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	ISO 15223-1, Clause 5.1.4	Use by date	Indicates the date after which the medical device is not to be used.
	ISO 15223-1, Clause 5.2.7	Not sterile	Indicates a medical device that has not been subjected to a sterilization process
	ISO 15223-1, Clause 5.7.7	Medical device	Indicates the item is a medical device
	ISO 15223-1, Clause 5.1.8	Importer	Indicates the entity importing the medical device into the locale
	ISO 15223-1, Clause 5.1.2	Authorized representative in the European Community	Indicates the authorized representative in the European Community.
	Medical Device Regulation 2017/745	CE marking	'CE marking' means a marking by which a manufacturer indicates that a device is in conformity with MDR 2017/745 and other applicable Union harmonisation legislation
	UK MDR 2002	UK Conformity Assessed Marking	Indicates that a device is in conformity with UK MDR 2002
	21 CFR 801.109	Prescription only	Requires a prescription in the United States.



Channel Medsystems, Inc.  
2919 7th Street  
Berkeley, CA 94710, USA

[www.cerene.com](http://www.cerene.com)  
[support@cerene.com](mailto:support@cerene.com)  
+1-510-338-9301



UK Responsible Person

MedEnvoy Global B.V.  
Prinses Margrietplantsoen 33  
Suite 123  
2595 AM The Hague  
The Netherlands

MedEnvoy UK Limited  
85, Great Portland Street, First Floor  
London, W1W 7LT  
United Kingdom



Kalms Consulting GmbH  
Rheinstr. 45-46  
12161 Berlin  
Germany  
cs@kalmsoperations.com  
+49 30 405 045 320



© 2023 Channel Medsystems, Inc. All rights reserved.  
*Cerene* is a registered trademark of Channel Medsystems.  
*Tyvek* is a registered trademark of E.I. du Pont de Nemours and Company or its affiliates.